## 2019-nCoV Ag Rapid Test Kit (Immunochromatography)

**Catalogue Number:**  
0555C2X001 (1 Test/Kit) 0555C2X005 (5 Tests/Kit) 0555C2X025 (25 Tests/Kit)

### INTENDED USE

This kit is used for the in vitro qualitative detection of 2019-nCoV antigen. It is an immunochromatography sandwich assay intended for the qualitative detection of 2019-nCoV nucleocapsid protein antigen in nasopharyngeal swab or nasal swab specimens. This kit is for professional use only.

### SUMMARY AND EXPLANATION

The novel coronaviruses belong to the  $\beta$  genus. COVID-19 is an acute respiratory infectious disease. People are generally susceptible. Currently, the patients infected by the novel coronavirus are the main source of infection; asymptomatic infected people can also be an infectious source. Based on the current epidemiological investigation, the incubation period is 1 to 14 days, mostly 3 to 7 days. The main manifestations include fever, fatigue and dry cough. Nasal congestion, runny nose, sore throat, myalgia and diarrhea are found in a few cases.

### PRINCIPLE OF THE TEST

This reagent uses double-antibody sandwich to legally detect the antigen of novel coronavirus (2019-nCoV) in nasopharyngeal swab or nasal swab samples. During detection, the gold labeled anti-2019-nCoV monoclonal antibody in the labeling pad binds to the 2019-nCoV antigen in the sample to form a complex, and the reaction complex moves forward along the nitrocellulose membrane under the action of chromatography, it is captured by the anti-2019-nCoV monoclonal antibody pre-coated in the detection zone (T) on the nitrocellulose membrane, and finally a red color reaction line is formed in the T zone. If the sample does not contain 2019-nCoV antigen, a red color reaction line cannot be formed in the (T) zone. Regardless of whether the sample to be tested contains 2019-nCoV antigen, a red reaction line will always form in the quality control area (C).

### MATERIALS AND COMPONENTS

#### Materials provided with the test kits

REF	0555C2X001	0555C2X005	0555C2X025
Component			
Test cassette	1	5	25
Swab	1	5	25
Tube with buffer	1	5	25

Instructions for use	1	1	1
Quick reference instructions	NA	1	1
Workstation	NA	NA	1

**Note:** The components in different batches of the kit cannot be mixed.

#### Materials required but not provided

1.Timer

### STORAGE AND STABILITY


- At 2°C - 30°C, the test device remains stable in the sealed pouch until the expiry date printed on the package. Do not store below 2°C and do not use expired product.
- The test cassette should be used within 15 minutes of removal from the foil pouch.
- The lot number and expiry date are printed on the labels.

### SAMPLE REQUIREMENTS

**Note:** The test cassette, the buffer tube and the specimen must be brought to room temperature (18°C - 26°C) before testing.


#### Nasopharyngeal swab:

- Insert the sterile swab into the place where the nasopharyngeal secretions are the most.
- Rotate the swab close to the inner wall of the nasal cavity 3 times.
- Remove the swab.

#### Nasal swab:

- Insert the entire absorbent tip of the swab (2-3 cm) into the left nostril and swipe firmly against the inside of the nostril more than 5 times in a circular motion.
- Take out the swab, insert it into the right nostril and firmly stroke the inside of the nostril more than 5 times in a circular motion.
- Take out the swab.


#### Note:

- False negative results may occur if the swab is not properly collected.
- The samples should be used as soon as possible after collected (within half an hour).
- Samples should not be inactivated.

### TEST PROCEDURE

#### Sample Processing:


- Open the drip cap of the buffer tube.
- Immerse the swab tip with the sample first into the buffer tube and rotate the swab tip 10 times in the buffer liquid while applying pressure with your fingers.
- Remove the swab while squeezing the sides of the tube to extract the liquid with the sample from the swab.
- Screw the drip cap tightly onto the buffer tube.
- Break off the tip of the drip cap to place 3 drops in the sample well of the test cassette.

### TEST PROCEDURE

Place the test cassette flat on the table and mark it.

- Add 3 drops or 80µL of the sample solution to the sample well and start the timer. It is recommended to use a pipette to collect buffer/sample to reduce any variation.
- As the test begins to work, you will see a colored line move across the result window in the center of the test cassette.
- Read the result at 15 minutes. Do not read the results after 20 minutes.

### INTERPRETATION OF TEST RESULTS

This test can only perform qualitative detection of 2019-nCoV antigen.

#### Positive Result:


If both C and T lines are visible at 15 minutes, the test result is positive and valid.

#### Negative Result:

If test area (T) has no colored line and the control area (C) displays a colored line, the result is negative and valid.

#### Invalid Result:

The test result is invalid if a colored line does not form in the control area (C). The sample must be re-tested, using a new test cassette.

### LIMITATIONS

- The result of the test should not be taken as a confirmed diagnosis, for clinical reference only. Judgement should be made along with RT-PCR results, clinical symptoms, epidemiological information and further clinical data.
- The contents of this kit are to be used for the qualitative detection of 2019-nCoV antigens from nasopharyngeal swab or nasal swab.
- Test performance depends on the amount of virus (antigen) in the sample and may or may not correlate with viral culture results performed on the same sample.
- The Sample buffer and test cassette must be equilibrated to room temperature (18°C~26°C) before used, otherwise the results may be incorrect
- A negative test result may occur if the level of antigen in a sample is below the detection limit of the test or if the sample was collected or transported improperly.
- Failure to follow the Test Procedure may adversely affect test performance and/or invalidate the test result.
- Reading time less than 15 minutes may lead a false negative result; Reading time more than 15 minutes may lead a false positive result.
- Positive test results do not rule out co-infections with other pathogens.
- Negative test results are not intended to rule in other viral or bacterial infections.
- Negative results should be treated as presumptive and confirmed with a molecular assay.
- Clinical performance was evaluated with frozen samples, and performance may be different with fresh samples.
- Users should test specimens as quickly as possible after specimen collection.
- If the sample volume is not enough, the chromatography cannot be carried out successfully. It is recommended to use a pipette to add samples.

### PERFORMANCE CHARACTERISTIC

#### 1. Clinical Verification

#### Nasopharyngeal swab performance:

The performance of 2019-nCoV Ag Rapid Test Kit (Immunochromatography) was compared with 356 RT-PCR nasopharyngeal swab samples, including 116 Positive samples with COVID-19 symptoms within seven days after onset of symptoms, 124 symptomatic Negative samples and 116 Asymptomatic Negative samples collected from donors.

2019-nCoV Ag Rapid Test Kit (Immunochromatography)	Comparative RT-PCR Test Result		
	Positive (+)	Negative (-)	Total
Detected Positive	113	1	114
Detected Negative	3	239	242
Total	116	240	356
Sensitivity	97.41%, 95% CI (92.67, 99.12)		
Specificity	99.58%, 95% CI (97.68, 99.93)		
Accuracy	98.88%, 95% CI (97.15, 99.56)		

The performance of 2019-nCoV Ag Rapid Test Kit (Immunochromatography) with positive results stratified by the comparative method cycle threshold (Ct) counts were collected and assessed to better understand the correlation of assay performance to the cycle threshold. As presented in the table below, the positive agreement of the 2019-nCoV Ag Rapid Test Kit (Immunochromatography) is higher with samples of a Ct count <25.

2019-nCoV Ag Rapid Test Kit (Immunochromatography)	Comparative RT-PCR Method (Positive by Ct Value)	
	Positive (Ct≤25)	Positive (Ct>25)
Detected Positive	85	28
Total	85	31
Positive agreement	100%	90.32%

**Nasal swab performance:**

The performance of 2019-nCoV Ag Rapid Test Kit (Immunochromatography) was compared with 350 RT-PCR nasopharyngeal swab samples, including 230 COVID-19 RT-PCR negative samples and 120 COVID-19 RT-PCR positive samples. In the evaluation, positive samples are collected both from symptomatic and asymptomatic groups. 104 samples out of 120 positive samples are collected from the patients that have one of the following symptoms within 7 days: headache, sore throat, difficult in breath, flu alike symptoms, back pain, body pain, etc. 210 of the 230 samples from the negative samples are collected from the asymptomatic PCR confirmed negative groups, 20 samples were collected from symptomatic patients.

2019-nCoV Ag Rapid Test Kit (Immunochromatography)	Comparative RT-PCR Test Result		
	Positive (+) (Ct<36)	Negative (-)	Total
Detected Positive	115	1	116
Detected Negative	5	229	234
Total	120	230	350
Sensitivity	95.83%, 95%CI (90.60%, 98.21%)		
Specificity	99.57%, 95%CI (97.58%, 99.92%)		
Accuracy	98.29%, 95%CI (96.31%, 99.21%)		

The performance of 2019-nCoV Ag Rapid Test Kit (Immunochromatography) with positive results stratified by the comparative method cycle threshold (Ct) counts were collected and assessed to better understand the correlation of assay performance to the cycle threshold. As presented in the table below, the positive agreement of the 2019-nCoV Ag Rapid Test Kit (Immunochromatography) is higher with samples of a Ct count <32.

Positive results broken down by Ct value:

2019-nCoV Ag Rapid Test Kit (Immunochromatography)	Comparative RT-PCR Method (Positive by Ct Value)
	Positive (Ct<32)
Detected Positive	115
Total	116
Positive agreement	99.14%

**2. Limit of Detection**

When the virus culture concentration was 100 TCID<sub>50</sub>/mL and above, the positive rate was greater than or equal to 95%. At virus culture concentration of 50 TCID<sub>50</sub>/mL and below, the positive rate is not higher than 95%, so the minimum detection limit of the 2019-nCoV Ag Rapid Test Kit (Immunochromatography) is 100 TCID<sub>50</sub>/mL.

**3. Cross-reactivity**

Cross-reactivity of the Kit was evaluated. The results showed no cross reactivity with the following specimen.

No.	Specimen type	Conc.
1	HCoV-HKU1	10 <sup>5</sup> TCID <sub>50</sub> /mL
2	Staphylococcus aureus	10 <sup>6</sup> CFU / mL
3	Streptococcus pyogenes	10 <sup>6</sup> CFU / mL
4	Measles virus	10 <sup>5</sup> TCID <sub>50</sub> /mL
5	Paramyxovirus parotitis	10 <sup>5</sup> TCID <sub>50</sub> /mL
6	Adenovirus 3	10 <sup>5</sup> TCID <sub>50</sub> /mL
7	Mycoplasma pneumoniae	10 <sup>6</sup> CFU / mL
8	Parainfluenza virus 2	10 <sup>5</sup> TCID <sub>50</sub> /mL
9	Human Metapneumovirus (hMPV)	10 <sup>5</sup> TCID <sub>50</sub> /mL
10	Human coronavirus OC43	10 <sup>5</sup> TCID <sub>50</sub> /mL
11	Human coronavirus 229E	10 <sup>5</sup> TCID <sub>50</sub> /mL
12	Human coronavirus NL63	10 <sup>4</sup> TCID <sub>50</sub> /mL
13	MERS-Coronavirus EMC/2012	10 <sup>4</sup> TCID <sub>50</sub> /mL
14	Bordetella parapertussia	10 <sup>6</sup> CFU / mL
15	Influenza B (Victoria strain)	10 <sup>5</sup> TCID <sub>50</sub> /mL
16	Influenza B (Y strain)	10 <sup>5</sup> TCID <sub>50</sub> /mL
17	Influenza A (H1N1 2009)	10 <sup>5</sup> TCID <sub>50</sub> /mL
18	Influenza A (H3N2)	10 <sup>5</sup> TCID <sub>50</sub> /mL
19	Avian influenza virus (H7N9)	10 <sup>5</sup> TCID <sub>50</sub> /mL
20	Avian influenza virus (H5N1)	10 <sup>5</sup> TCID <sub>50</sub> /mL

21	Epstein-Barr virus	10 <sup>5</sup> TCID <sub>50</sub> /mL
22	Enterovirus CA16	10 <sup>5</sup> TCID <sub>50</sub> /mL
23	Rhinovirus	10 <sup>5</sup> TCID <sub>50</sub> /mL
24	Respiratory syncytial virus	10 <sup>5</sup> TCID <sub>50</sub> /mL
25	Streptococcus pneumoniae	10 <sup>6</sup> CFU / mL
26	Candida albicans	10 <sup>6</sup> CFU / mL
27	Chlamydia pneumoniae	10 <sup>6</sup> CFU / mL
28	Bordetella pertussis	10 <sup>6</sup> CFU / mL
29	Pneumocystis jirovecii	10 <sup>6</sup> CFU / mL
30	Mycobacterium tuberculosis	10 <sup>6</sup> CFU / mL
31	Legionella pneumophila	10 <sup>6</sup> CFU / mL

**4. Interference Substances**

The test results do not be interfered with the substance at the following concentration:

No.	Interference substances	Conc.
1	Whole Blood	4%
2	Ibuprofen	1mg / mL
3	Tetracycline	3µg / mL
4	Chloramphenicol	3µg / mL
5	Erythromycin	3µg / mL
6	Tobramycin	5%
7	Throat spray (Menthol)	15%
8	Mupirocin	10mg/mL
9	Throat lozenge (Menthol)	1.5mg/mL
10	Tamiflu (Osetamivir)	5mg/mL
11	Naphthoxoline hydrochloride nasal drops	15%
12	Mucin	0.50%

13	Fisher's Friend	1.5mg/mL
14	Compound Benzocaine Gel	1.5mg/mL
15	Cromoglycate	15%
16	Sinex (Phenylephrine Hydrochloride)	15%
17	Afrin (Oxymetazoline)	15%
18	Fluticasone propionate spray	15%

**5. Precision**

- Test replicates of negative and positive by using the reference materials of enterprises. The negative agreement and the positive agreement were 100%.
- Test three different lots kits including positive and negative reference materials of enterprises. The negative results and the positive results were 100%.

**6. Hook Effect**

There was no Hook effect detected when the concentration of inactivated virus stock solution raised up to 4.0×10<sup>5</sup> TCID<sub>50</sub>/ml.

**PRECAUTIONS**

- For in vitro diagnostic use.
- Do not use the kit contents beyond the expiration date printed on the outside of the box.
- Use appropriate precautions in the collection, handling, storage, and disposal of patient samples and used kit contents.
- Use of Nitrile, Latex (or equivalent) gloves is recommended when handling patient samples.
- Do not reuse the used test cassette, tubes or swabs.
- The user should never open the foil pouch of the test cassette exposing it to the ambient environment until the test cassette is ready for immediate use.
- Discard and do not use any damaged or dropped test cassette or material.
- The sample buffer contains a salt solution (saline). If the solution contacts the skin or eye, flush with copious amounts of water.
- Inadequate or inappropriate sample collection, storage, and transport may yield false test results.
- Sample collection and handling procedures require specific training and guidance.
- Use the appropriate Fixed Volume Pipette in accordance with test procedures.
- To obtain accurate results, do not use visually bloody or overly viscous samples.
- To obtain accurate results, an opened and exposed test cassette should not be used.
- Testing should be performed in an area with adequate ventilation.
- Wear suitable protective clothing, gloves, and eye/face protection when handling the contents of this kit.
- Wash hands thoroughly after handling.

**KEY TO SYMBOLS USED**

Consult instructions for use

Date of manufacturer

Store at 2°C~30°C

Do not reuse

Use-by date

Keep away from Sunlight

Manufacturer

Contains sufficient for <n> tests

Batch code

Catalogue number

In vitro diagnostic medical device

Keep dry

Authorized representative in the European Community

**Guangzhou Decheng Biotechnology Co., LTD**

Room 218, Building 2, No.68, Nanxiang Road, Science City, Huangpu District, 510000, Guangzhou P.R.China  
Tel.: +86-020-82557192  
service@dochebio.com  
www.dochebio.com

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**Caretechion GmbH**